



Unomedical Pty Ltd
11-17 Wilmette Place,
Mona Vale, NSW, Australia 2103
T (02) 9997 8033
F (02) 9999 4031

New Zealand
Customer Service: 0800 161 777
Product Support: 0800 445 400
www.unomedical.com

Anaesthesia



Local response, global expertise

Unomedical leadership

As you would expect from Unomedical, their Anaesthesia product range provides patient-considerate, purpose-designed circuits and ancillaries carefully developed from only highest quality to use. High performance at realistic prices.

Anaesthesia Products

- Anaesthetic Circuits
- Obstetric Circuits
- Accessories
- Filters
- ET/TT Tubes




Anaesthetic Circuits

Anaesthetic Circuits are manufactured in Australia. Unomedical is able to custom manufacture circuits to meet customer requirements

Product Features

- High quality corrugated tubing with good drapability and flexibility
- Kink resistant
- Connectors (22mm) conform to AS/NZS 2496 and provide secure leak-resistant attachment
- Available in a range of lengths and configurations
- Neoprene (non latex) breathing bags
- Single use design reduces risk of contamination
- Custom options available



Corrugated Circuit Safety Design

No blockages will occur if the circuit is accidentally twisted or kinked.

REF.	Description	Order Qty
1675-00	Swivel Wye without Sampling Port	50
1750-00	30 Metre Corrugated Tube Cuffed every 150mm - Blue	1
1791-15	1.5 Metre Circuit with Gas Sampling Port & 2 L Neoprene Reservoir Bag with 1.5 Metre Tube and extra 1.2 Metre tube (4 pieces)	20
1791-30	3.0 Metre Circuit with Gas Sampling Port & 2 L Neoprene Reservoir Bag with 3.0 Metre Tube and extra 1.2 Metre tube (4 pieces)	20
1795-15	1.5 Metre Circuit with Gas Sampling Port & 2 L Neoprene Reservoir Bag with 1.5 Metre Tube	20
1795-18	1.8 Metre Circuit with Gas Sampling Port & 2 L Neoprene Reservoir Bag with 1.5 Metre Tube	20
1795-27	2.7 Metre Circuit with Gas Sampling Port & 2 L Neoprene Reservoir Bag with 1.5 Metre Tube	15
1800-18	1.8 Metre Circuit, Wye Swivel without Gas Sampling Port	25
1781-00	1.8 Metre Circuit, Wye Swivel with Gas Sampling Port	25
1796-18	1.8 Metre Circuit with Gas Sampling Port and 2L Neoprene Reservoir Bag with 1.5 Metre Tube and Soft Cuff Connector	20

Obstetric Circuits



2.7 METRE OBSTETRIC CIRCUIT WITH MOUTH PIECE

REF.	Carton Size	Weight	Qty
1788-27	548 x 382 x 567	5.1kg	25



2.7 METRE OBSTETRIC CIRCUIT WITH MOUTH PIECE AND FILTER

REF.	Carton Size	Weight	Qty
1790-27	548 x 382 x 567	5.5kg	20

Accessories

Manufactured from safe, stable non-toxic material. Designed to extend range of intended use.



REF.	Description	Carton Size	Weight	Qty
1670-00	22mm T Piece	200 x 165 x 250	1kg	50
1671-00	Mouth Piece	200 x 165 x 250	1kg	100
1674-00	Swivel Wye with Sampling Port	200 x 165 x 250	1.1kg	50

Gas Sampling Line

NON STERILE



REF.	Description	Qty
656.10.300	With 3 Metre Male Luer Lock Connector	100

Filters

SlimLine Breathing Filter Range



6200/100 HMEF

6100/01 Multi Pleat

6000/100 Multi Pleat Plus HMEF



6200/01

9070/01 Paediatric

9066/01 Paediatric HMEF

REF.	Filter Media Type	Humidification (mg/L H ₂ O)	Deadspace (ml)	Resistance to flow	Weight	Gas port Efficiency	Bacterial Filtration Efficiency	Viral Filtration Volume	Recommended tidal	Qty/ Box
SLIMLINE MULTI PLEAT PLUS WITH LUER PORT										
6000/01	Pleated/ Hydrophobic	23	80	3.1cm H ₂ O 60 LPM	33	yes	99.99999	99.9999	>200ml	
SLIMLINE MULTI PLEAT PLUS HMEF WITH LUER PORT										
6000/100	Pleated/ Hydrophobic	28	80	3.9cm H ₂ O 60 LPM	33	yes	99.99999	99.9999	>200ml	50
6000/100-S*										
SLIMLINE MULTI PLEAT WITH LUER PORT										
6100/01	Pleated/ Hydrophobic	19.7	42	4.2cm H ₂ O 60 LPM	29	yes	99.9999	99.9999	>150ml	50
6100/01-S*										
SLIMLINE WITH LUER LOCK										
6200/01	Electrostatic	x	42	3.1cm H ₂ O 60 LPM	22	yes	99.9997	99.992	>150ml	50
6200/01-S*										
SLIMLINE HMEF WITH LUER LOCK										
6200/100	Electrostatic	33.1	42	3.7cm H ₂ O 60 LPM	22	yes	99.9997	99.992	>150ml	50
6200/100-S*										
SLIMLINE PAED HMEF WITH LUER LOCK										
9070/01-S*	Electrostatic	31	27	1.4cm H ₂ O 30 LPM	16	yes	99.9997	99.997	>100ml	50
SLIMLINE NEONATAL HMEF WITH LUER LOCK										
9080/100-S*	Electrostatic	30	9	0.3kpa @ 5 LPM	10	yes	99.999	99.998	>50ml	50
SLIMLINE NEONATAL WITH LUER LOCK										
9080/01-S*	Electrostatic	x	9	0.2kpa @ 5 LPM	9	yes	99.999	99.998	>50ml	50

Note: * = Sterile

ET/TT Tubes


Regular



Regular – High Volume Low Pressure Cuff

REGULAR MURPHY

REF.	Size I.D. (mm)	Qty / Box
MM6110050	5.0	10
MM6110055	5.5	10
MM6110060	6.0	10
MM6110065	6.5	10
MM6110070	7.0	10
MM6110075	7.5	10
MM6110080	8.0	10
MM6110085	8.5	10
MM6110090	9.0	10



Regular – Standard Volume Low Pressure Cuff


REGULAR MURPHY

REF.	Size I.D. (mm)	Qty / Box
MM61150050	5.0	10
MM61150055	5.5	10
MM61150060	6.0	10
MM61150065	6.5	10
MM61150070	7.0	10
MM61150075	7.5	10
MM61150080	8.5	10
MM61150085	8.5	10
MM61150090	9.0	10
MM61150095	9.5	10


Regular – Plain (Without Cuff)

REGULAR MURPHY

REF.	Size I.D. (mm)	Qty / Box
MM61130020	2.0	10
MM61130025	2.5	10
MM61130030	3.0	10
MM61130035	3.5	10
MM61130040	4.0	10
MM61130045	4.5	10
MM61130050	5.0	10
MM61130055	5.5	10
MM61130060	6.0	10
MM61130065	6.5	10



Intubating Stylet



REF.	Size FR.	O.D. mm	Qty / Box
MM64200006	6	2.0	10
MM64200010	10	3.3	10
MM64200014	14	4.9	10

ET/TT Tubes

Curved

Oral Curved – Plain (Without Cuff)



ORAL CURVED MURPHY


REF.	Size I.D. (mm)	Qty / Box
MM61233030	3.0	10
MM61233035	3.5	10
MM61233040	4.0	10
MM61233045	4.5	10
MM61233050	5.0	10
MM61233055	5.0	10
MM61233060	6.0	10
MM61233065	6.5	10
MM61233070	7.0	10

Oral Curved – High Volume Low Pressure Cuff



ORAL CURVED MURPHY

REF.	Size I.D. (mm)	Qty / Box
MM61213050	5.0	10
MM61213055	5.5	10
MM61213060	6.0	10
MM61213065	6.5	10
MM61213070	7.0	10
MM61213075	7.5	10
MM61213080	8.0	10
MM61213085	8.5	10
MM61213090	9.0	10



Nasal Curved – HighVolume Low Pressure Cuff

NASAL CURVED MURPHY

REF.	Size I.D. (mm)	Qty / Box
MM61313060	6.0	10
MM61313065	6.5	10
MM61313070	7.0	10
MM61313075	7.5	10
MM61313080	8.0	10

Nasal Curved – Plain (Without Cuff)



NASAL CURVED MURPHY

REF.	Size I.D. (mm)	Qty / Box
MM61333030	3.0	10
MM61333035	3.5	10
MM61333040	4.0	10
MM61333045	4.5	10
MM61333050	5.0	10
MM61333055	5.5	10
MM61333060	6.0	10
MM61333065	6.5	10
MM61333070	7.0	10